

# Light Curing Unit

## — Description

Huge Light Curing Unit is designed for the polymerization of light curing individual tray and other light curing materials. It is featured with high polymerizing capacity and low maintenance. The shaped light curing material can be fast cured in 3 minutes in it. And it is quite easy to operate with just a few steps in need.



## — Advantages

- Easy operation
- High polymerization capacity
- All voltage available: 110V and 220V

## — Instructions



Lay the Light Curing Unit to be horizontal and connect it to power supply



Put the shaped light curing material into the chamber for curing



Switch on the Unit, choose the 180s mode, press the button "Start" and wait a moment for curing



Easy operation · High accuracy · High plasticity · Pleasant smell

## Light Curing Tray / Unit

— An Ideal Custom Tray Material

Preformed to save your precious time

GLOBAL DENTAL SUPPLIER



Shandong Huge Dental Material Corporation

Add / No.68 Shanghai Road, Donggang District, Rizhao City, Shandong Province, 276800, P.R. China.

Tel / +86 (633) 2277268 marketing@hugedental.com www.hugedental.com

Facebook Huge Dental

Instagram Huge Dental

Youtube Huge Dental

# Light Curing Tray

## – Description

Huge Light Curing Tray is a preformed tray indicated for making individual tray and can also work as temporary base plate in dental labs. As light curing material, it can be fast cured in 3 minutes in a unit. The product is characterized by easy operation, high precision, low deformation rate, high plasticity and pleasant mint flavour.



## – Application

- Making individual tray



Types Features	Light Curing Tray	Self-cured resin Tray	Metal Tray
Operation	preformed → ready for use and time-saving	need mixing and kneading first → complicated and time-consuming	-
Plasticity	soft but with certain strength → easy to trim and shape when opening windows	too soft to work with in dough stage → difficult to trim and shape when opening windows	hard and rigid metal material → difficult to open window
Accuracy	high hardness after curing → low deformation rate ensuring accuracy	existence of contraction → unstable and inaccurate	specific and uniform size → fail to match different oral forms and inaccurate
Smell	mint flavour → pleasant smell	residual stimulative monomer → pungent smell	-

- Making temporary base plate

Types Features	Light Curing Base Plate	Wax base Plate
Operation	preformed → ready for use and time-saving	need heating first → time-consuming
Accuracy	high hardness after curing → low deformation rate ensuring accuracy	vulnerable to temperature → unstable and inaccurate



## – Main Advantages

- Preformed → easy handling
- 3-minute fast curing → quick curing
- High hardness and fit different oral forms → accurate
- Soft → high plasticity
- Mint flavour → pleasant smell
- Various choices → pink or blue, thicker or thinner, all available

## – Specification

Standard package



50pcs/box

Sample package



2pcs/bag

## – Storage

Store this product in the ventilation. Avoid flame, high temperature and direct sunlight.

Storage temperature: 4°C ~ 25°C

Validity: 3 years