

PMMA Discs and Blocks

Excellent biocompatibility, shade stability and esthetics.

Product					
	Monolayer PMMA	Multilayer PMMA	Full Denture PMMA	Pink PMMA	Clear PMMA
Shades	16A-D, A0 B0, BL, BL1-4	16A-D, 2M2, BL, BL1-4	Denture Base: 2S, 2ST, 3S, 3ST, Pink Brown Denture Arch: 16A-D, 2M2, BL, BL1-4	2S, 2ST, 3S, 3ST, Pink Brown	Clear, Pink, Orange, Red
Flexural Strength	>120Mpa	>120Mpa	>120Mpa	>120Mpa	>120Mpa
Modulus of Elasticity	>2200Mpa	>2200Mpa	>2200Mpa	>2200Mpa	>2200Mpa
Water Absorption	<20µg/mm ³	<20µg/mm ³	<20µg/mm ³	<20µg/mm ³	<20µg/mm ³
Water Solubility	<0.8µg/mm ³	<0.8µg/mm ³	<0.8µg/mm ³	<0.8µg/mm ³	<0.8µg/mm ³



Shandong Huge Dental Material Corporation

Add / No.68 Shanghai Road, Donggang District, Rizhao City, Shandong Province,
276800, P.R. China.

Tel / +86 (633) 2277268 marketing@hugedental.com www.hugedental.com

[Facebook](#) Huge Dental

[Instagram](#) Huge Dental

[Youtube](#) Huge Dental